

Bath Safety

Shower Chairs – Stationary

- Your delivery technician will assemble the unit.
- Place the shower chair completely in tub or shower stall. Check for stability – all feet flat on floor of bathing area.
- Raise seat height to desired seat height.
 - In a shower: match to a chair seat height that is comfortable for you.
 - In a tub: raise seat height to top of tub side or higher to make transfers easier.

Tub Transfer Benches

- These models are intended for those who will transfer from a wheelchair to the bathtub. They can be used in shower stalls, but you will need to modify the cover (shower curtain) over the opening to contain water from the shower spray.
 - Place the bench half in the tub (2 legs in). The backrest should be on the edge away from the faucets. (User seat should face faucets). The feet fitted with the suction cups should be inside the tub.
 - Adjust the seat height by depressing the snap button lock and moving the leg extensions up or down the legs of the bench. Note: Depending on your bathroom's design, the legs outside the tub may need to be set at a different level than the two legs inside the tub.
- ** Make sure the seat is level and stable before using.